



**HUMAN CARE**  
Your Life. Your Way.

**The new range of  
true floor level beds**

# *Introducing the new true floor level beds*

*The cost of fall related injuries is huge for both health care systems and the society. Fall injuries are without doubt the biggest problem in hospital care today. Serious fall injuries lead to extended hospital stay and increased cost for hospitals or care providers. These injuries also have physical and psychological effects on patients, from severe pain and loss of confidence after a bed fall to permanent loss of mobility.*

Ten centimeters (four inches). That is the difference between falling – or rolling – out of bed! So small a distance, yet so large a difference between extended hospital stays and cost – or their avoidance. Human Cares new range of true floor level beds lowering to 10 cm (4”), can to reduce fall injuries and provide the opportunity for excellent care.

Patients who are at risk of falling are not protected from the impact of a fall from a usual low bed that only descends to around 20 cm (8”). That distance to the floor – although appearing small and harmless – can still lead to serious injuries.

The use of floor level beds also greatly enhances the dignity of care. Functioning at a floor level, fall injuries can be prevented without the use of restraints. Side rails can deprive patients of their dignity and autonomy and also impose a risk of injuries due to entrapment.

For an optimum working environment for caregivers the Human Care range of floor-level beds can be raised to any ergonomically suitable height up to 799 mm (31,4”) from the floor, helping to minimize workplace back injury risks.

*” 10 cm (4 inches) makes the difference between falling or rolling out of bed – These 10 cm avoid extended hospital stays and direct costs. ”*

# range of

## FloorLine-*i* Plus – a complete care management solution

The aesthetically appealing FloorLine-*i* Plus is not just a bed, it is a complete care management solution. With its extensive range of features such as trendelenburg, cardiac chair position and electronic CPR function, the FloorLine-*i* Plus offers the best possible care. With its open bus digital control box it is also future proofed for adding even more accessories.

The FloorLine-*i* Plus is also the perfect solution for preventing fall injuries. Its minimum height from the floor is just 99 mm [3,9"]. At the same time no compromises are made for excellent ergonomic working height reaching 799 mm [31,4"].



## FloorLine-*i* – a sound purchase

The FloorLine-*i* is a sound purchase. It offers a variety of functionality such as trendelenburg and cardiac chair position providing excellent care.

The FloorLine-*i* is an optimal floor-level-bed solution for use both in hospitals and long term care facilities. Its minimum height from the floor is just 99 mm [3,9"] and at the same time no compromises are made for caregivers ergonomics.



## FloorLine LTC – your choice for long term care

The FloorLine LTC with its home environment design is the optimal choice for long term care facilities. Its one-piece design makes it easy to fold and store away efficiently.

The FloorLine LTC is the perfect solution for preventing fall injuries. Its minimum height from the floor is just 95 mm [3,7"] providing a safe and considerate way of taking care of people who may be at risk of falling out of their bed.



Technical Specifications	FloorLine- <i>i</i> Plus	FloorLine- <i>i</i>	FloorLine LTC
Part no.	Europe: 90500 Australia 1: 90501 Australia 2: 90502 USA: 90503 United Kingdom: 90504	Europe: 90300 Australia : 90301 USA: 90302 United Kingdom: 90303	Australia : 90608 Australia with decorative panel : 90609 USA : 90604 USA with decorative panel: 90605 USA class 1 : 90606 USA class 1, with decorative panel: 90607
Platform position	99 - 799 mm / 3.9 - 31.4"	99 - 799 mm / 3.9 - 31.4"	95 - 700 mm / 3.7 - 27.5"
Sleeping dimension	2000 (2175) x 900 mm / 78.7 (85.6) x 35.4 "	2000 (2175) x 900 mm / 78.7 (85.6) x 35.4 "	2000 (2175) x 900 mm / 78.7 (85.5) x 35.5 "
Bed dimension	2325 (2500) x 925 mm / 91.5 (98.4) x 36.4 "	2325 (2500) x 925 mm / 91.5 (98.4) x 36.4 "	2150 (2325) x 905 mm / 84.5 (91.5) x 35.5 "
User Weight	185 kg / 408 lbs	185 kg / 408 lbs	150 kg / 330 lbs
Safe working load	250 kg / 551 lbs	250 kg / 551 lbs	200 kg / 440 lbs



# HUMAN CARE

Your Life. Your Way.

## SWEDEN

Årstaängsvägen 21C  
117 43 Stockholm  
Phone:- +46 8 665 35 00  
Fax:- +46 8 665 35 10  
[info.se@humancaregroup.com](mailto:info.se@humancaregroup.com)

## AUSTRALIA

Unit 9, 271 Wells Road  
Chelsea Heights, Victoria 3196  
Phone:- +61 3 8773 1111  
Fax:- +61 3 9773 4008  
[info.au@humancaregroup.com](mailto:info.au@humancaregroup.com)

## UNITED STATES

4210 S. Industrial Drive, Suite 160  
Austin, TX 78744  
Phone:- 512.476.7199  
Fax:- 512.476.7190  
[info.us@humancaregroup.com](mailto:info.us@humancaregroup.com)

## CANADA

10-155 Colonnade Road  
Ottawa, ON K2E 7K1  
Phone:- 613.723.6734  
Fax:- 613.723.1058  
[info.ca@humancaregroup.com](mailto:info.ca@humancaregroup.com)