




HUMAN CARE

Your Life. Your Way.



Convertible Chairs



"I'm no different from you. I might be older, less agile. My legs are slow but my mind is quick. I have a career. Or maybe I've had one. Either way I'm not slowing down. I can keep up. But some things don't come easy anymore. Small things, and some big ones. Appearance matters, I'm still vain. I'm still happy. I know I need some help, even more than just a little. I don't mind as long as it works.

But hey, it is what it is."

Your Life. Your Way.

Our focus is simple. We work for our users. Our products serve to be a natural part of life to anyone who might need them. Meeting that need is what drives and motivates us.

We are highly experienced in medical care and mobility aids. With our dedicated and knowledgeable sales team, we always guarantee high quality service and reliable support.

Welcome to Human Care.



Convertible Chairs

Our Convertible Chair series offers a multitude of capabilities within a single product. Designed to allow a single care-giver to easily and safely transport and reposition a patient, the chairs can perform the functions of stretcher, mobile chair, and custom seating.

H-250

- Provides safe transfer to and from bed
- Limitless number of positions
- Single caregiver transfers made simple



The H-250 Convertible Chair is a unique device that easily converts from a stretcher to an upright mobile chair in seconds. Its Tilt-in-Space capability allows for increased blood circulation helping to prevent skin breakdown and other complications of immobility.

With it's 113 kg/250 lbs capacity, the H-250 provides a no-lift solution for moving, repositioning and transferring patients in the safest manner possible.

The H-250 Convertible Chair offers the ability to move between a bed and a chair with ease and enables a single caregiver to easily transfer a patient.

*Safe, comfortable
and convenient.*

H-250

Part no.	Seat Width	Overall width	Headrest adjustable height	Maximum user weight	Accommodates varying patient heights
B1002-BL	431 mm / 17"	610 mm / 24"	178 mm / 7"	113 kg / 250 lbs	1270 - 1980 mm / 5' 0" - 6' 6"

I-400

- Sturdy but lightweight
- Unlimited repositioning capabilities
- Easy to turn and maneuver



Our I-Series of convertible chairs has an unlimited number of reclined and tilt positioning. The ability to frequently reposition a user significantly reduces the chance of pressure ulcers and increases the users tolerance for long-term seating.

A sturdy frame can be equipped with a patient transfer system to enable a single caregiver to perform safe lateral transfer to and from bed.

The I-Series Chairs have three different weight capacities to suit all needs; 181 kg/400 lbs, 317 kg /700 lbs and 454 kg/1000 lbs.

Perfect for repositioning and transfer.



I-700

- Tilt-in-space function
- Helps to reduce risk of injuries
- Fully adjustable head support

I-1000

- Ideal for bariatric patients
- Full range of optional accessories
- Minimizes pressure ulcers and sores



I-400, I-700 and I-1000

Product	Part no.	Seat Width	Overall width	Overall length	Maximum user weight	Accommodates varying patient heights
I-400	B1041-TS-GEL	546 mm / 21.5"	749 mm / 29.5"	1880 mm / 74"	181 kg / 400 lbs	Manual Tilt and Recline
I-700	B1071-TS-GEL	686 mm / 27"	889 mm / 35"	1880 mm / 74"	318 kg / 700 lbs	Manual Tilt and Recline
I-1000	B1083-TS-GEL	838 mm / 33"	986 mm / 38.6"	1880 mm / 74"	454 kg / 1000 lbs	Electric tilt and recline



HUMAN CARE
Your Life. Your Way.

AUSTRALIA

Unit 9, 271 Wells Road
Chelsea Heights, Victoria 3196
Tel:- +61 3 8773 1111
Fax:- +61 3 9773 4008

Email: info.au@humancaregroup.com

CANADA

10-155 Colonnade Road
Ottawa, ON K2E 7K1
Tel:- 613.723.6734
Fax:- 613.723.1058

Email: info.ca@humancaregroup.com

SWEDEN

Årstaängsvägen 21C
117 43 Stockholm
Tel:- +46 8 665 35 00
Fax:- +46 8 665 35 10

Email: info.se@humancaregroup.com

UNITED STATES

4210 S. Industrial Drive, Suite 160
Austin, TX 78744
Tel:- 512.476.7199
Fax:- 512.476.7190

Email: info.us@humancaregroup.com