

# QUATTRO SW



## Stackable Wheelchair



QUATTRO SW provides the perfect modularity for transfer with comfort, speed, easy handling, mobility and space saving features. With its fantastic design QUATTRO SW is convenient to use in hospitals, airports and shopping centers.

### MESPA HEALTHCARE PRODUCTS INC.

2. organize sanayi bolgesi 83216 nolu cadde no:20 Sehitkamil-Gaziantep / TURKIYE

Tel: +90 342 337 5310

Fax:+90 342 337 53 13

e-mail: [export@mespa.com.tr](mailto:export@mespa.com.tr)

web: [www.mespa.com.tr](http://www.mespa.com.tr)



- To enhance safety, QUATTRO SW has an integrated brake system that is automatically activated as the brake handle is released
- Stainless steel structure gives QUATTRO SW a longer service life along with easy clean property
- ABS footrest section is foldable upward when idle and stacking
- Horizontally stacking 10 chairs into 3, 10 meters, QUATTRO SW brings a new dimension to storage
- Backrest, seat surface and armrest made of ABS to provide seat comfort together with hygiene
- Height adjustable stainless steel I.V. pole
- 360° swivel small front castors enhance maneuverability
- Urine bag holder

#### Optional:

- *Coin lock system*
- *Safety belt*



TECHNICAL DETAILS	
Depth:	900mm
Width:	641mm
Height:	1000mm
Front Castors Diameter:	125mm
Rear Castors Diameter:	200mm
Backrest Dimensions:	455x300mm
Seat Padding Dimensions:	460x405mm
Footrest Dimensions:	260x250mm
Safe Working Load Capacity:	150kg
Weight:	14,5kg

### MESPA HEALTHCARE PRODUCTS INC.

2. organize sanayi bolgesi 83216 nolu cadde no:20 Sehitkamil-Gaziantep / TURKIYE

Tel: +90 342 337 5310

Fax: +90 342 337 53 13

e-mail: [export@mespa.com.tr](mailto:export@mespa.com.tr)

web: [www.mespa.com.tr](http://www.mespa.com.tr)

